

Power cleanse. Product boost. One device.



- Dual-action power for your best skin experience.
- Skinvigate™ Duo features soft, silicone bristles and vibration technology that is clinically shown to effectively dislodge dirt, oil and makeup.
- Take your skin care routine to the next level with thermal massage which is known to help temporarily dilate pores and support better product absorption.

Did you know?

Your face has 43 muscles and the ability to make more than 20 facial expressions. That's at least 43 reasons to pair your favorite Mary Kay® facial oil, moisturizer, boosters or serum with the thermal massage function to give your face the special treatment it deserves.*

Customized Experience: Features five intensities to allow you to control cleansing and thermal massage with the touch of a button.

Skinvigorator™ Duo is a 2-in-1 facial device that elevates your routine with deep cleansing and thermal massage to enhance your next skin care steps. Designed for smart multitasking, Skinvigorator™ Duo features dual-head vibration technology: On one side, soft silicone bristles pulse to cleanse the skin, while the warming massage side helps boost product absorption. With the press of a button, you can switch from the cleansing side to the thermal massage side — no replaceable device heads to swap out! The beauty of double-duty skin care starts here, because the best skin care routines deserve the best tools.

Two Possibilities, One Device:

- **Elevate the Everyday Cleanse:** Did you know that hand cleansing alone can leave pore-clogging grime behind? Clinical data shows that vibration technology effectively dislodges dirt, oil and makeup and improves circulation. Skinvigorator™ Duo features the gentle pulsing of soft, silicone bristles that act like dozens of tiny fingers hard at work to help cleanse your skin.
- **Warmth That Works Wonders:** Take your skin care routine to the next level with thermal massage. Heat is known to help temporarily dilate pores and support better absorption of skin care products. With Skinvigorator™ Duo, applying your facial oil, moisturizer, boosters or serum* becomes a pampering experience, enhanced by soft pulses and a warming function.

Application Instructions:

Cleansing function:

Remove eye makeup first. Wet face and facial cleansing brush head with water, then apply your favorite Mary Kay® cleanser directly to your face. Turn on device and select the desired vibration intensity. Gently cleanse using swiping or massaging motions with the device. Avoid eye area and hairline. After use, hold the power button down to turn off the device. Rinse face thoroughly. Rinse facial device head with water, and allow to air dry or gently dry with a soft towel. If needed, use a mild soap to clean Facial Cleansing Brush Head. Rinse and air dry. Suitable for daily cleansing, morning and evening. DO NOT submerge in water. Thoroughly dry the device before placing it back in the stand. Please read all instructions carefully before use. DO NOT use in the shower. *Tip – select a lower vibration intensity in the morning for a gentle cleanse, and a higher vibration intensity in the evening to remove makeup.*

Thermal massage function:

Apply facial oil, moisturizer, boosters or serum* to a dry face. Turn on the device and select the desired vibration intensity. Press the warm function button to turn on the thermal massage feature. Carefully apply the heating surface to the face. Use gentle swiping or massaging motions. Avoid the eye area and hairline. After use, hold the power button down to turn off device. Rinse facial device head with water, and allow to air dry or gently dry with a soft towel. Suitable for daily use. DO NOT submerge in water. Thoroughly dry the device before placing it back in the stand. Please read all instructions carefully before use. DO NOT use in the shower.

Device features:

- Water-resistant.
- Silicone bristles offer a hygienic alternative to traditional nylon brushes — rinsing clean in seconds with no need for replacement heads.
- Suitable for daily use on all skin types, including sensitive skin.
- Rechargeable battery with USB cord means no searching for replacement batteries.
- White light indicates device is charging. Green light indicates device is fully charged.
- Charging time is approximately three hours.

WARNING: To reduce the risk of electrocution: Keep the power plug (not included), USB cord and charging port dry. Lay the device flat to dry before using the stand. Do not plug/unplug the device with wet hands. When charging, do not place or store the device near water (e.g., bathtub, shower, sink, toilet). Do not use the device while charging. Discontinue use if the product is damaged or not working properly.

CAUTION: For external use only. Use only as directed. Do not use in bathtub or shower. Do not use damaged parts or parts not designed for the Skinvigorator™ Duo Facial Device. Keep away from children and pets. This is not a toy. To prevent suffocation, keep packaging materials, including plastic bags, out of reach of babies and children. Do not use on irritated or damaged skin. If redness or irritation develops, stop use. If the redness or irritation persists, contact your healthcare provider. Store in a cool, dry place away from direct sunlight.

*The device was not designed to be paired with acne, retinol or sunscreen products.